

Accessories for Western Blotting

Solutions for your Western Blotting workflow

Azure offers a variety of accessories for electrophoresis and Western blotting.

Blotting Paper for Western Blotting

Pre-cut blotting papers save time and reduce waste for the easiest Western blotting workflow. Forget about measuring and cutting blotting paper to the size of your mini-gels. Use these pre-cut papers to assemble your transfer “sandwiches”. Available in two different sizes.



Blot Development Folders

A wrinkle-free alternative to plastic wrap, these transparent plastic folders hold your blot during handling and imaging and help maximize the signal longevity of Radiance and Radiance Plus substrates.

Background Quenching Sheets

Use these quenching sheets under your blot when imaging chemiluminescent or fluorescent Westerns. The sheets reduce background due to contamination or auto-fluorescence, improving the quality of your image. Available in two different sizes.



Incubation Trays

These specially designed incubation trays are ideal for staining and washing electrophoresis gels and membranes. Available in a range of sizes in both clear and opaque.



Advantages

- Lidded – to protect experiments from dust or debris
- Non-binding – surfaces treated with silicone to prevent antibody loss
- Convenient – multiple sizes available for different applications
- Opaque or clear – Opaque trays are ideal for dealing with light sensitive fluorescent blots

Ordering Information

Part Number	Name	Description	Size
Blotting Paper for Western Blotting			
AC2110	Blotting Paper for Western Blotting	Pre-cut blotting paper for Western blot transfer – 7x9 cm vertical gels	60 sheets
AC2111	Blotting Paper for Western Blotting	Pre-cut blotting paper for western blotting transfer – 8x10 cm vertical gels	60 sheets
Blot Development Folders			
AC2126	Blot Development Folders	Transparent plastic sheets to hold blots for imaging – alternative to plastic wraps	100 sheets
Background Quenching Sheets			
AC2144	Quenching Sheets	Background quenching sheets, for use with fluorescent or chemiluminescent samples	10 sheets
AC2147	Quenching Sheets, Large	Background quenching sheets for use with fluorescent or chemiluminescent samples – 16x22 cm	10 sheets
Incubation Trays			
AC2120	Incubation Tray – 7x5 cm (opaque black top and bottom)	For staining and washing light sensitive gels and membranes	5 trays
AC2121	Incubation Tray – 9x6 cm (opaque black top and bottom)	For staining and washing light sensitive gels and membranes	4 trays
AC2122	Incubation Tray – 11x9 cm (opaque black top and bottom)	For staining and washing light sensitive gels and membranes	3 trays
AC2123	Incubation Tray – 15.5x10.5 cm (opaque black and bottom)	For staining and washing light sensitive gels and membranes	2 trays
AC2150	Incubation Tray – 7x5 cm (clear)	For staining and washing gels and membranes	5 trays
AC2151	Incubation Tray – 9x6 cm (clear)	For staining and washing gels and membranes	4 trays
AC2152	Incubation Tray – 11x9 cm (clear)	For staining and washing gels and membranes	3 trays
AC2153	Incubation Tray – 15.5x10.5 cm (clear)	For staining and washing gels and membranes	2 trays



www.azurebiosystems.com • info@azurebiosystems.com

Copyright © 2016 Azure Biosystems. All rights reserved. The Azure Biosystems logo and Azure™ are trademarks of the Company. All other trademarks, service marks and trade names appearing in this brochure are the property of their respective owners.

DC0039-001