

Document: D-54073  
Version: D-54073-001  
Release date: Jul 1, 2014  
Update date: Jul 1, 2014  
Printed on: Jul 1, 2014

Applies to part number: AC2127

# MATERIAL SAFETY DATA SHEET

## Section 1 Product and company identification

Product name: **Azure Transfer Buffer**  
Supplier: Azure Biosystems, Inc.  
Address: 6761 Sierra Court, Suite F  
Dublin, CA 94568  
(925) 307-7127

Recommended Use: For research use only.

Restrictions on Use: Not for clinical use. Not for internal use in animals or humans. Not for diagnostic use. Not for household or any other unintended use.

## Section 2 Hazards identification

### EMERGENCY OVERVIEW:

Harmful if swallowed. May cause eye and skin irritation. May cause respiratory tract irritation. May be harmful if absorbed through skin. Avoid breathing vapors, or dusts. Use with adequate ventilation. Avoid contact with eyes, skin, and clothes. Wash thoroughly after handling. Keep container closed.

HMIS RATING	
HEALTH	2
FLAMMABILITY	0
REACTIVITY	0

NFPA RATING	
HEALTH	2
FLAMMABILITY	0
REACTIVITY	0

For additional information on toxicity, please refer to Section 11.

## Section 3 Composition/Information on Ingredients

Ingredient	CAS #	%
Tris-(hydroxymethyl)-aminomethane	77-86-1	3
Tris-(hydroxymethyl)-aminomethane acetate	6850-28-8	9
Disodium ethylenediaminetetraacetate dihydrate	6381-92-6	1
Sodium dodecyl sulfate	151-21-3	<1

## Section 4 First Aid Measures

ORAL EXPOSURE	If swallowed, wash out mouth with water provided person is conscious. Call a physician.
INHALATION EXPOSURE	If inhaled, remove to fresh air. If breathing becomes difficult, call a physician.
DERMAL EXPOSURE	In case of contact, immediately wash skin with soap and copious amounts of water.
EYE EXPOSURE	In case of contact with eyes, flush with copious amounts of water for at least 15 minutes. Assure adequate flushing by separating the eyelids with fingers. Call a physician.

## Section 5 - Fire Fighting Measures

FLASH POINT	N/A
AUTOIGNITION TEMP	N/A
FLAMMABILITY	N/A
EXTINGUISHING MEDIA	
Suitable:	Water spray. Carbon dioxide, dry chemical powder, or alcohol-resistant foam.
Not suitable:	Not known
FIREFIGHTING	
Protective Equipment:	Wear self-contained breathing apparatus and protective clothing to prevent contact with skin and eyes.

## **Section 6 Accidental Release Measures**

METHODS FOR CLEANING UP: Absorb on sand or vermiculite and place in closed containers for disposal. Ventilate area and wash spill site after material pickup is complete.

## **Section 7 - Handling and Storage**

### HANDLING

User Exposure: Avoid inhalation. Avoid contact with eyes, skin, and clothing. Avoid prolonged or repeated exposure.

### STORAGE

Suitable: Keep tightly closed. Store at room temperature.

## **Section 8 - Exposure Controls / PPE**

### PERSONAL PROTECTIVE EQUIPMENT

Respiratory: Respiratory protection not required. If desired, use multi-purpose respirator.

Hand: Protective gloves.

Eye: Chemical safety goggles.

GENERAL HYGIENE MEASURES: Wash thoroughly after handling.

## **Section 9 - Physical and chemical properties**

APPEARANCE PHYSICAL STATE: Colorless, clear liquid. Water-based, pH 8

## **Section 10 Stability and Reactivity**

STABILITY: Stable under recommended storage conditions

Materials to Avoid: Strong oxidizing agents, acids, bases.

### HAZARDOUS DECOMPOSITION PRODUCTS

Hazardous Decomposition Products: Carbon oxides, nitrogen oxides.

### HAZARDOUS POLYMERIZATION

Hazardous Polymerization: Will not occur.

## **Section 11 - Toxicological Information**

### ROUTE OF EXPOSURE

Skin Contact: Causes skin irritation.

Skin Absorption: May be harmful if absorbed through the skin.

Eye Contact: Causes eye irritation.

Inhalation: May be harmful if inhaled. Causes respiratory tract irritation.

Ingestion: Harmful if swallowed.

### SIGNS AND SYMPTOMS OF EXPOSURE

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

## **Section 12 Ecological Information**

No known significant effects or critical hazards.

## **Section 13 Disposal Considerations**

### APPROPRIATE METHOD OF DISPOSAL OF SUBSTANCE OR PREPARATION:

Observe all federal, state, and local environmental regulations.

Within the present knowledge of the supplier, this product is not regarded as hazardous waste.

## **Section 14 Transport Information**

### DOT

Proper Shipping Name: None

Non-Hazardous for Transport: This substance is considered to be non-hazardous for transport.

### IATA

Non-Hazardous for Air Transport: Non-hazardous for air.

## **Section 15 Regulatory Information**

**Section 16 Other information**

Date of last revision: May 9, 2014.

For R&D use only. Not for drug, household or other uses.

To the best of our knowledge, the information contained herein is accurate. However, neither the above-named supplier, nor any of its subsidiaries, assumes any liability whatsoever for the accuracy or completeness of the information contained herein. Final determination of suitability of any material is the sole responsibility of the user. All materials may present unknown hazards and should be used with caution. Although certain hazards are described herein, we cannot guarantee that these are the only hazards that exist.